

55\*115

## Intelligent Wireless Dental Camera

### User Manual



#### Please read carefully this user manual before turn on device

Welcome use the intelligent visual appartus dental camera. In order to protect your security and profits, please check product list and carefully ready user manual before use device.

#### 1 Precautions

To avoid the risk of malfunction, electric shock, injury, fire, and damage to equipment or property, be sure to observe the following safety precautions.


1. This product has a built-in battery. Do not throw this product into a fire or expose it to the sun, otherwise the battery may explode.
2. Do not use this product while charging.
3. Do not modify or repair this product by yourself.
4. Do not put the charging cable of this product in a humid environment or water
5. Do not place this product where children and infants can reach.

6. Do not drop or impact/shock the product.
7. Do not use this product while bathing or soaking in hot springs.
8. Before and after use, wipe and disinfect the entrance part with an alcohol cotton pad.


#### 2 Product parameter

Product name: <b>Intelligent Visual Dental Camera</b>	Product Model: <b>T6</b>
Charging mode: <b>Type-C</b>	Input parameters: <b>5V-500mA</b>
Camera Resolution: <b>1MP</b>	Working current: <b>200mA</b>
Working power: <b>0.88W</b>	Operating Voltage: <b>3.56-4.2V</b>
Battery Capacity <b>550mAh</b>	Standby time: <b>2H</b>

#### 3 Button/Port Explanation



⏻ ON/OFF	Long press Button
📷 Take Photos	Click Button
📹 Video Record	Click Button
💡 Light Adjustabl	Click Button to Switch 4 Gears



#### 4 Usage

- 1.APP downloading  
For Android:  
1) Search"Timesiso"on Google Play Store  
2) Scan the QR code below to download and install "Timesiso"  
For IOS:  
Search"Timesiso"on Apple APP Store to install the APP



2. Long press the ON/OFF key, and the WIFI indicator flashes on the upper left of the screen.
3. Open the APP,the indicate window pop out,click OK to jump to Mobile Settings - WLAN
4. Connect the WiFi named "Timesiso-T6XXXX" with the phone, and the password is 12345678. After successful connection, the indicator light will turn on.
5. Open the software, and you can see the images and videos taken by the camera on your phone.
6. Stick the camera in the month and adjust the Angle and brightness.

7. Press the corresponding button to complete the photographing and recording operation.
8. After finishing the photographing, you can open the album to view the photographing results.
9. After use, please follow the following steps:
  - 1) Disconnect WIFI;
  - 2) Long press the ON/OFF key to shut down;
  - 3) Clean the device with an alcohol cotton pad and keep it in a cool and ventilated place.

#### 5 Malfunction solutions

- 1.No response when turn on/off,the camera doesn't light up  
Solutions: Please try again after full power
- 2.APP flash back  
Solutions: Clean up the running background of mobile phone , restart the APP
- 3.The red lamp inside the dental camera blinking  
Solutions: Indicate the device is in low power, please charge after use.

- 4.Failed to display the real-time image properly after the device successfully connected to APP  
Solutions: Disconnect the WIFI connection between the phone and the endoscope, restart the device or clean the background, open the APP again and reconnect to WIFI.

- 5.WiFi connection is unstable or no images appear  
Solutions: Restart the device and reconnect to the APP  
\*If the above solutions still fail to solve the problem, please contact customer service, we will further provide you with practical solutions.

#### 6 Charging

When the power is low, the red light will flash; you need to use the power adapter to charge, using the specified 5V/1A.  
When the battery is charging, the indicator light is red.  
After the charging, the indicator light turns into green (charging for 2 hours).  
When the battery is fully charged, the product can be used for about 1.5 hours.

#### 7 After sales policy

- 1.Within half a year from the date of purchase,if the device has quality problems

due to human factors,can contact customer service ask for maintenance,Please fill out the registration form truthfully,and return the defective products back to us.

- 2.Within half a year from the date of purchase,if the product cannot be used due to its own problems,Can contact customer service ask for exchange,According to the process of exchange and return,We'll send the new device out within 7 days after we receive the defective device.

- 3.Warranty NOT covered for the following conditions,Our Company or our authorized service center will collect reasonable fee for maintenance from you:

- 1)Exceed the period of validity of "three guarantees".
- 2)No valid invoice or altered warranty certificate
- 3)Causing losses due to not use or maintain according to the instructions.
- 4)Damage the devices due to destroy or repair without the authorization of our company.
- 5)Soaking,breaking,or the printed circuit board burned out.
- 6)Caused damage to the surface of the device when using,beside the shell damage due to the structural or material factors.
- 7)Damage caused by force majeure