
尺寸：5.5cm\*11.5cm

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CN

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
EN

IMPORTANT HINT

1. Please fully charge the device before using it for the first time, and charge it for at least 1 hour before it can be used normally. Please do not charge directly through the PC terminal. Please choose a 5V 1A adapter.

2. If you use a microscope to observe cells or microorganisms, the material used for observation with the microscope must be thin and transparent, so you need to make your own specimen slides for observation.

3. This equipment cannot accurately read the magnification of the microscope. It is a microscope that combines digital and optical magnification. The specific magnification effect is subject to the actual picture taken.

4. When connecting to the app, please click allow to read the location information and storage space. Your private information will not be read.

5. The microscope cannot be used on both the mobile phone and the PC at the same time.

6. Do not use mobile phone network and home wifi during use.

7. When the device is connected to the computer to view images, the camera button of the device cannot be used.

8. The best focal length of the microscope is 0-70mm. You need to adjust the focal length through the focal length adjustment roller to achieve the clearest state.


9. WiFi connection is only applicable to your mobile phone and tablet, PC computer cannot be used, if you want to use it on PC computer, please connect via a USB data cable, and please download the correct computer software to use.

10. Please close the useless APP in your phone to ensure that the microscope runs smoothly and will not get stuck or crash.

11. Do not touch the lens and other optical parts with your hands, otherwise the image will be blurred and the image quality will be affected.

12. Do not disassemble the product or component structure yourself to avoid abnormal use of the equipment.

13. If there are static black spots when the screen is displayed, please wipe the lens plate with a cleaning cloth until the black spots disappear, as shown in the figure below.

PRODUCT INTRODUCTION

Thank you for purchasing our WiFi electron microscope, which can be easily used in different fields, including:

1. Textile industry for textile inspection

2. Printing inspection

3. Industrial inspection: PCB, precision machinery

4. Educational purpose

5. Hair check

6. Skin examination

7. Microbial observation

8. Jewelry and coins (collection) check

9. Visual aids

10. Other

This is a portable WiFi electron microscope with WiFi hotspot. It not only supports connection to iOS/Android system mobile phones and tablets, but also supports USB connection to the computer. The larger the screen, the better the display effect and the clearer the image quality.

The device can take pictures, video and file storage through mobile APP and computer software.

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


(二) Bracket fixed

After assembling the lens, place the host into the bracket positioning hole and tighten the bracket fixing nut clockwise. Place the stand on a stable surface (an unstable surface will cause the picture to shake).

(三) Image viewing

Mobile users


1. APP download and installation

a. iOS:

Search for "WiFi Check" in the App Store to download and install, and then use the device.


b. Android:


Search for "WiFi Check" in Google Play to download and install.

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Quick focus of lens


Quick focus of lens

Note: When using the No. 1 lens to observe the object under test, the stage can only move up and down within the No. 1 adjustment range marked by the bracket. When observing the object under test with No. 2, the stage can only move up and down within the No. 2 adjustment range marked by the bracket, and no imaging is allowed beyond the corresponding range.

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Through the fourth step of operation, we can determine the approximate object distance, and then adjust the actual multiple required according to the actual situation. The bracket adjustment and the host adjustment are corresponding. If it is in a clear state, the bracket stage is adjusted downwards. The corresponding focus adjustment roller needs to be rotated to the right, and the multiple is reduced. The bracket stage is adjusted upwards, and the corresponding focal length adjustment roller needs to be rotated to the left, and the multiplier increases.

7. Introduction and use of mobile APP interface

Take a photo: Click to take a photo.


Video: Click to record video/end video recording.

Save files: Click to view or delete the photos or videos that have been taken.

Rotation: Click to rotate the observation object 90°.

Mirror: Click the image to flip left and right.

Resolution: Click to switch between different resolutions.

——PC user

\*Note: When using a computer, the camera button of the device cannot be used.

(MAC users)

1. For Macbook (Mac OS X 11.0+)

Please directly use the Photo Booth or Quick Time Player bundled with Macbooks.

2. In the "Applications" directory of the Finder window, find an application named "Photo Booth".

3. Connect the device

Connect to the computer via a USB cable, then long press to turn on the device and run "Photo Booth". Select the camera "GENERAL-UVC", as shown in the figure below. Ready to use (when connected to the computer, the green indicator light will flash slowly).

(Windows users)

1. For win 10, you can directly insert the USB port and search for CAMERA in WINDOWS to find it and click.

For win 7/8, downloading "Amcap" will solve this problem.


Note: Please disable the default laptop camera in Windows! You must change the privacy settings of the camera, and you need to make changes to allow access.

2. Or download and install the software "Smart Camera" from the following link: [www.inskam.com/download/camera.zip](http://www.inskam.com/download/camera.zip)

3. Connect the device

Connect to the computer via a USB cable, and then long press to turn on the device.

When you plug it into the computer for the first time, it will automatically install the USB driver. Please wait for the computer to successfully display the installation. If the installation fails, please plug it in again or change another USB port.

PRODUCT PARAMETER

Pixel	5,0M pixels
Magnification	10X-2000X
WiFi connection distance	5m (open environment)
Working height	0-70mm
Focus mode	Manual focus
Knob focus range	2mm-55mm
Dynamic frame number	20-200/f/s
USB interface type	USB2.0
Data format	JPG/MP4
Light source	10 bright LED lights +1 super bright LED light (brightness adjustable)
Operating system	Android 8.0, iOS 9.0 or above Windows XP, Vista, Win7/8/10, Mac 11.0 or above
Power supply	18650 lithium battery 3.7V/1800mAh Working time: 5 hours; Charging time: 3 hours

TROUBLESHOOTING

Problem	Solution
Device cannot be turned on	1. Is the battery charged? 2. Is power connected 3. Is it charged for more than 1 hour when there is no electricity
No image on phone	1. Check the WiFi settings on the phone to confirm whether the device is successfully connected. 2. Close the APP and reopen it. 3. Turn off the device and turn it on again. 4. Check if the phone allows WiFi Check to connect to wireless. 5. Check whether the mobile phone has closed the data network.
The image is not clear	1. Confirm whether the measured object is correctly positioned directly under the lens. 2. Whether the focal length has been adjusted to the best position 3. Whether the light is used correctly
Dark spots appear on the image	1. Clean the lens lens carefully with a soft cloth 2. Whether there is dirt on the surface of the measured object 3. Whether the stage is dirty
The microscope is connected to the PC, After successfully downloading and opening the software, there is no image	1. Reconnect the microscope and PC 2. Connect to the USB port on the back of the host computer 3. Whether the battery has low power alarm 4. Is the data cable used correctly?
Device crashes	1. Press the reset hole of the host with a fine needle to restart the device
The battery indicator does not light red when charging	1. Check whether the charging head is normal 2. Check whether the cable is normal 3. Check if the battery is properly assembled
Screen image has black edges	1. Is the lens fixed in the middle of the bracket 2. Whether the object being observed is placed upright
Camera feedback (IOS device)	1. Go to the "Settings" page and find "WiFi Check" in the application section at the bottom to allow "Photos" to be turned on.