


Intelligent Visual Ultrasonic Dental Scaler



User manual
Model: T11 Pro
EN

1 Product description

Please carefully read user manual before use for correct operation.
This product is a new type of visual ultrasonic dental scaler. The product is characterized by connecting mobile phones and tablets and other devices to intuitively understand the condition of your teeth through the App. It uses ultrasonic high-frequency resonance technology to remove stains on the teeth and achieve the effect of whitening and protecting teeth.



Button function definition:
Power button:
Long press to turn on/off, short press to switch gears
Short press to switch gears:
Low-end/mid-end/high-end
LED:
Power on: Green light flashes
Connect to wifi: the green light is always on.
LED:
Low battery: red light flashes
Charging: the red light is always on
Full: light off

EN-2

2 Safety Precautions

To avoid the risk of malfunction, electric shock, injury, fire and damage to equipment or property, please follow the following safety precautions

1. This product has a built-in battery. Do not throw this product into a fire or expose it to the sun, otherwise it may cause the battery to explode.
2. Do not use this product while charging.
3. Do not use this product when the charging port is not properly covered.
4. Do not modify or repair this product by yourself.
5. Do not place the charging money of this product in a humid environment or in water.
6. Do not place this product where children and infants can reach.
7. Do not drop or impact/shock the product.
8. Do not use this product while bathing or bathing in hot springs.
9. Those who suffer from severe periodontitis, whose teeth are being treated, or who are worried about symptoms of discomfort in the oral cavity, please consult a doctor before use.
10. Do not use paint thinner, benzene or other chemicals to clean this product.
11. When using this product for the first time, there may be slight bleeding, because the root of the tooth is exposed to this kind of irritation for the first time, which is normal. If there is severe bleeding, please stop using this product immediately. Please also consult a doctor.

EN-3

3 Steps

Step one:
Scan the QR code with your mobile phone and download the APP (Android4.4, iOS9.0 Or the above system can support the use)



EN-4

Step two:

- 1 Open the APP, click on the upper right corner + Add device
- 2 Bind the visual scaler device; press and hold the product for three seconds to turn on, and check the above steps to complete. Click to start connecting.



EN-5

Step three:

Click to start connection and jump to mobile phone settings-wireless LAN.
Connect to WIFI: SUNUO-T11-XXXX Password: 12345678



EN-6

3 Product stall introduction

Caueat:


This product belongs to the working principle of ultrasonic. After turning on, it will not produce over-frequency vibration when it touches the gums or soft objects, and it will have obvious working performance when it touches teeth or hard objects to produce over-frequency resonance!

According to your needs, select the mode gear you need.

Milder Mode: Soft power for sensitive teeth, small tartar stains

Clean Mode: Normal power for stubbon tartar such as smoke stains, coffee stains, tea stains, etc.

Powerful Mode: Strong power for deep cleaning such as formed dental calculus, hard tartar stains.



EN-7

Frequently Asked Questions Answers

1. Is gum bleeding normal during use?

Answer: If people with periodontitis or gingivitis use it or it is easy to cause gum bleeding, including brushing teeth or biting hard food, it will also bleed. If the gums are bleeding when using this product, just stop using it and rinse with cold water.

2. Is it normal for the tip of the cleaning tooth to become hot during use?

Answer: When using, electrical energy is converted into kinetic energy, and it is normal for the temperature of the scaler to rise slightly.

3. Is it normal that the product sometimes has sound?

Answer: The sound at work is caused by the resonance of the working tip and the teeth, which is a normal phenomenon.

4. Is it normal for the product to suddenly stop working during ultrasonic work?

Answer: The product is equipped with an automatic shutdown function, which will stop the work of the transducer after the transducer has been working for 10 minutes, and you need to switch the gear again to start the working state.

EN-8

Troubleshooting

1. There is no response when booting, and the camera does not light up

Fully charge before use.

2. APP crashes

Clean up the background of the mobile phone and restart the APP.

3. The red light of the scaler flashes

Prompt this product is low battery, please use it after charging.

4. The product fails to display real-time images normally after successfully connecting to the Sunuo App

Disconnect the WIFI connection between the mobile phone and the scaler, restart the machine or clean the background, reopen the App, and reconnect to the WIFI.

5. The WiFi connection is unstable or no image appears.

Restart the device and reconnect to the APP

* If the problem is still not solved according to the above operations, please contact customer service, we will further provide you with practical solutions.

EN-9

Product parameter

Product name Intelligent visual ultrasonic dental scaler	Product number T11 Pro
Network standard IEEE 802.11b/g/n	Working frequency 2.4GHz
Image transfer rate 30fps	Image Sensor CMOS
Working temperature -10~50 °C	Battery 1200mAh
Battery life -About 90 minutes	Charging time About 1.5h -2h
Input parameters 5v == 0.5A	Operating Voltage 5V
Working power 1.68W	Pixel 5MP
Best focus 0~2CM	

Made in China
CE FC RoHS
EN-10