

**SUNUO**

## Visual Wifi Ultrasonic Tooth Cleaner

Best selling Sunuo® T12 pro Visual wifi ultrasonic tooth cleaner with camera, you can see intraoral situation when use the teeth cleaning tools, Visual wifi ultrasonic tooth cleaner made in China, Built-in 2 million pixels high-definition camera, LED lights and WIFI technology, so that you can not see the condition of your teeth clearly when using the device, and check your teeth status easily. WiFi connected electric calculus remover makes it easier to use the product, download the app in the app store of your iOS and Android devices. You can download the app on Google Play. Once you turn on the power of the product, you can connect it to your phone via WiFi.

## Product Description

### FEATURES:

- 1.5 Megapixels HD endoscopic camera: high-definition endoscope camera can accurately locate the dirt location.
- 2.Two kinds of professional cleaning heads to meet your different cleaning needs, perfect cleaning of tooth dirt in different positions.
- 3.Very fast WIFI image chip. Watching while cleaning, the image transmission is smooth and does not freeze.
- 4.Three gears professional cleaning modes, according to the dirt of different difficulty for precise cleaning.
- 5.IPX7 waterproof, the whole machine is IPX7 waterproof and can be washed directly.
- 6.The product have smart sensor, When the cleaning head touches the teeth, it will automatically start cleaning. When it touches the gums, the cleaning head stops working to protect your teeth.



Product parameters	
Product name	Intelligent VISUAL ultrasonic dental scaler
Product model	T12 Pro
Network stanard	IEEE 802.11b/g/n
Working frequency	2.4Ghz
Image transfer rate	30 fps
Image Sensor	CMOS
Best focus	0-2cm

Camera pixel	5 megapixels
Product color	White/ Blue/ Black
Working temperature	-10℃~50℃
Battery life	About 90 minutes
Charging time	About 1.5h-2h
Input parameters	DC 5V/ 0.5 A
Working power	1.8W
Battery capacity	1200mAh
Waterproof rate	IPX7
Charge cable	Type-C

**SUNUO® T12 Pro Visual Wifi Ultrasonic Tooth Cleaner**



**Package Included**

- 1x Ultrasonic scaler
- 2x Clean the needle
- 2x Protective cover
- 1x Charge cable
- 1x Wrench
- 1x User manual

1x Packing box



### How to connect the phone?

- 1: Download SUNUO APP in Google Store
- 2: Open the APP, click on the upper right corner "+" Add device
- 3: Bind the visual scaler device: press and hold the product for three seconds to turn on, and check the above steps to complete. Click to start connecting.
- 4: Click to start connection and jump to mobile phone settings-wirelessLAN.  
Connect to WiFi: SUNUO-T12 pro-XXXX Password: 12345678

