

**Disposable Medical Cap**

As the professional manufacture, we would like to provide you high quality Fealter® Disposable Medical Cap. And we will offer you the best after-sale service and timely delivery. Disposable medical cap: The blue material is made of PP+PE non-woven fabric. Breaking strength $\geq 30\text{N}$, water resistance $\geq 1.67\text{kPa}$, The material has high water resistance, high breaking strength, no peculiar smell and so on, and has the function of general isolation.

You are welcomed to come to our factory to buy the latest selling, low price, and high-quality Fealter® Disposable Medical Cap. We look forward to cooperating with you.

Product name: Fealter® Disposable Medical Cap

Intended use: General isolation of medical personnel at work

Specification of disposable medical cap: One size

Specification	Cap circumference length	Cap height
one size	62 cm	14 cm
Deviation	$\pm 2\text{cm}$	$\pm 2\text{cm}$

Storage: Fealter® Disposable medical cap should be stored in a dry, ventilated and non-corrosive environment, away from fire sources and inflammables.

Precautions:

1. Please strictly follow the usage and usage restrictions of this product when using it. Users should judge whether this product is suitable for specific occasions. If the usage exceeds the usage and usage restrictions of this product, we do not make any guarantee or assume any responsibility related to this.

2. Please avoid direct sunlight.

3. This product melts at 135°C . Please pay attention to keep away from open flames and heat sources when using it.

4. Not washing, washing will affect the protective effect.

5. Do not iron, machine wash or dry clean.

6. Please dispose it properly according to relevant national regulations on pollutant disposal after use.

Product validity: Storage period is 3 years.