



Non-Contact Infrared Thermometer Forehead

The Meiyunsheng® is a professional manufactures of Non-Contact Infrared Thermometer Forehead. Our professional expertise in manufacturing Non-Contact Infrared Thermometer Forehead has been honed over the past 10+ years. This Non-Contact Infrared Thermometer forehead is designed for measuring the temperature of human body especial forehead without contacting for infants/adults, and it can keep the record of temperatures as well.

Product Description

Non-Contact Infrared Thermometer forehead

You are welcomed to come to our factory to buy the latest selling, and high-quality Non-Contact Infrared Thermometer forehead. We look forward to cooperating with you.

1. Introduction of the Non-Contact Infrared Thermometer forehead

The Meiyunsheng® Non-Contact Infrared Thermometer forehead is with high-precision infrared sensor for fast and accurate measurement, and large LCD screen, bright backlight, clear display, for easy reading.

2. Parameter (Specification) of the Non-Contact Infrared Thermometer forehead

Effective Distance: 5-15cm

Accuracy: ±0.2°C

Memory Data: 32 Sets Reading Recall

Operating Temperature: 10 to 40° C (50 to 104° F)

Storage Condition: (-25°--55°) °C RH≤85% Power Supply: D.C.3V2AA Alkaline Battery

Size: 156 x 95 x 42.5mm

Net Weight: 123g

3. Feature And Application of the Non-Contact Infrared Thermometer forehead



This Non-Contact Infrared Thermometer forehead is convenient for household or public use, ${}^{\circ}C/{}^{\circ}F$ selectable, and alarm in case of abnormal temperatures, auto data record and auto power off.

4. Product Details of the Non-Contact Infrared Thermometer forehead

LCD display of three backlight colors, buzzer sound selectable, data record of last 32 measurements. The default alarm point is $36.5^{\circ}\text{C}/97.7^{\circ}\text{F}$ and $38^{\circ}\text{C}/100.4^{\circ}\text{F}$.

5. Product Qualification of the Non-Contact Infrared Thermometer forehead

The Meiyunsheng® Non-Contact Infrared Thermometer forehead is certified with CE, FCC, FDA certifications etc.

6. Deliver, Shipping and Serving of the Non-Contact Infrared Thermometer forehead

We usually keep a certain stock of thisNon-Contact Infrared Thermometer forehead, for orders under 10 pcs, delivery with in 10 days, for that over 20k pcs, delivery within 20 days, and it could be shipped as per your requirement.