

Extended Version EmeTerm Fashion Blue Replacement Strap



Emeterm is a healthcare solution trading and manufacturer focusing on wearable and treatable electronic wrist device for anti-nausea and vomiting and extending to medical consumables. We consist of leaders in the field of medical research and product development who are devoted to provide effective and convenient therapy options to users

around the world now. We have over 10 years experience in medical business, approved with FDA, UKCA and ISO etc. EmeTerm Blue Replacement Strap - Extended Version exported to over 20 authorized distributors worldwide, covering markets in more than 30 countries or regions including the United States, Canada, the United Kingdom, Australia, and Japan. We are looking for more business partner and cooperation around the world .

Extended Version EmeTerm Fashion Blue Replacement Strap

Emeterm is a healthcare solution trading and manufacturer focusing on wearable and treatable electronic wrist device for anti-nausea and vomiting and extending to medical consumables. We consist of leaders in the field of medical research and product development who are devoted to provide effective and convenient therapy options to users around the world now. We have over 10 years experience in medical business, approved with FDA, UKCA and ISO etc. EmeTerm Blue Replacement Strap - Extended Version exported to over 20 authorized distributors worldwide, covering markets in more than 30 countries or regions including the United States, Canada, the United Kingdom, Australia, and Japan. We are looking for more business partner and cooperation around the world .

WAT Medical is a Canadian company that focuses on developing innovative and personalized medical technologies. The team consists of international medical experts, winning multiple national competitions and funding since 2016. The company has reached an estimated market value of 7 million USD. Since the global market expansion in 2018, the sales revenue has reached over \$1.4 million USD. The innovative products of WAT Medical consist of properties of both medical and consumer technologies.



EmeTerm Blue Replacement Strap - Extended Version introduction

EmeTerm Blue Replacement Strap - Extended Version is a safety and effective anti-nausea wristband that prevents and relieves nausea and vomiting induced by it can be used for various conditions of dizziness with zero side-effects. We provide a one-year warranty for all products. It is covered with waterproof technology and 30 hours battery time for outdoor activities. Electric Plates come out as a non-invasive neuronal electrical stimulation that goes through the brain to prevent brain transmitting nausea signals to the stomach, then reduces the nauseous sensation for the user. Immediately effective for the treatment of nausea and vomiting caused by motion sickness (while in a car, a boat, or a plane), PONV (post-operative nausea and vomiting) and morning sickness (during pregnancy) with zero side-effects. Our product is easy to install and charge. We have

helped a lot people suffering from motion sickness. They felt absolutely no nausea and enjoy life now.



Product Specifications:

Brand	EmeTerm
color	Black, Blue
Principle	TENS Targeted Neuromodulation
Indications	Effective relief of all forms of nausea& vomiting including morning sickness, sea sickness, car sickness, VR sickness, air sickness, vertigo.
Age Range	The EmeTerm device is indicated for the treatment of nausea and vomiting in patients 6 years of age or older.
Certificate	FDA, CE, TGA, HC, ISO
Award	IF Design 2017

Notice	Indications vary due to regulations. Please check with customer service and local suppliers.
--------	---

Success Rate:

US Version EmeTerm 1 Anti-nausea Wristband With an 85% success rate according to our clinical trials, the EmeTerm Anti-nausea Wristband is the FDA-cleared, all-natural wearable solution for nausea and vomiting.

Patented Conduction Enhancement Electrodes

Increased contact area

Better conduction

No gel required

Rechargeable Magnetic Charging

Equipped with a magnetic USB charging cable that can be quickly charged anytime and anywhere conveniently.

Only one charge is required for short trips.

HOW TO USE US Version EmeTerm 1 Anti-nausea Wristband?

STEP 1

Position: 2-3cm away from the transverse wrist line, with a width of about 2 fingers

STEP 2

Wear the device, Ensure skin is sufficiently in contact with electrodes

STEP 3

Turn on Home key to start the device, then adjust level until numb your palm.

Principle of US Version EmeTerm 1 Anti-nausea Wristband

Nausea and vomiting are caused by transmission of a nausea signal to the stomach via the human brain, which causes the stomach to contract irregularly, resulting in vomiting.

By releasing a particular low-frequency pulse, EmeTerm can adjust the vagus nerve signals traveling to and from the stomach, and can prevent or postpone it from receiving nausea signals. This then slows down the irregular gastric contraction, and reduces the nauseous sensation for the users.

1. TENS consists of small unit, attached to the wrist via electrode pads, delivers signal through the median nerve at the location on the underside of the wrist.
2. These signals travel through the body's nervous system to the part of the brain which controls nausea, which causes the stomach to contract irregularly, resulting in vomiting.
3. By releasing a low-frequency pulse, it can adjust the vagus nerve signal traveling to and from the stomach, and can prevent or postpone the brain from transmitting nausea signals to the stomach, then reduces the nauseous sensation for the users.