

















User Manual




Model No: CL-H001


Parts List:

<p>A</p>  <p>x1</p>	<p>B</p>  <p>x1</p>		
<p>C</p>  <p>x1</p>	<p>D</p>  <p>x1</p>	<p>E</p>  <p>x2</p>	<p>F</p>  <p>x1</p>
<p>G</p>  <p>x1</p>	<p>H</p>  <p>x1</p>	<p>I</p>  <p>5# x1</p>	<p>J</p>  <p>M4*8 x3</p>
<p>K</p>  <p>M4*6 x20</p>	<p>L</p>  <p>x1</p>	<p>M</p>  <p>x1</p>	<p>N</p>  <p>x1</p> <p>O</p>  <p>x4</p>

User Guide:




Mode: 1





Mode: 2

Mode:1




Mode: 2



Thickness of door 40MM



Mode: 2





Mode: 2



Mode: 2


3



Open

Health Mode

Parameter Display



Parameter setting



Lose fat

Lose weight mode



Toner skin

Toner skin mode



Fitness

Fitness mode



Nerve

Neural recovery mode



Sleep

Sleep mode



joint

Joint therapy mode



Wound

Trauma recovery mode




Animal

Animal therapy mode

Health Mode

Parameter Display



_____ Tap to set the Red Light, (+) increase and (-) to decrease.



_____ Tap to set the NIR, (+) increase and (-) to decrease.



_____ Tap to set the Pulse, (+) increase and (-) to decrease.



_____ Tap to set the TIME, (+) increase and (-) to decrease.



_____ Function increases.



_____ Function decreases

Hz

_____ Pulse

%

_____ Percentage of the set value



_____ Run/pause/press 5 seconds for the switch

Voice instruction content

Wake up the machine :

Hello Cutie

Mode operation instructions:

Lose fat please

Fitness please

Sleep please

Wound please

Tener skin please

Nerve please

Joint please

Animalplease

The parameters set for machine to operate:

Operation

Running machine to stop instruction:

Stop

Volume control instructions:

Volume up

Volume down

Voice function switch instructions:

Turn off voice function

Turn on voice function

Run time settings instructions:

Ten minutes

Twenty minutes

Thirty minutes

Forty minutes

Fifty minutes

Sixty minutes

Red light setting brightness instructions:

Red light zero

Red light ten percent

Red light twenty percent

Red light thirtypercent

Red light forty percent

Red light fifty percent

Red light sixty percent

Red light seventy percent

Red light eighty percent

Red light ninety percent

Red light one hundred percent

Infrared light setting brightness:

Infrared light zero

Infrared light ten percent

Infrared light twenty percent

Infrared light thirty percent

Infrared light forty percent

Infrared light fifty percent

Infrared light sixty percent

Infrared light seventy percent

Infrared light eighty percent

Infrared light ninety percent

Infrared light one hundred percent

Possible problems:

1.No voice response:

Check and see if the device is on mute,unmute it by command „Turn on voice function”

Check and see if the machine is running, it needs to be waked up by „Hello Cutie” and command to STOP”

Try again with more accurate pronunciation/speak more slowly.

Restart



2.The light is not on:please check whether the power supply is connected.

3.The red light doesn't hurt the eyes,but we suggest you wear an eye mask because it is bright.

Contact us

Website: www.errayhealing.com

E-mail: info@errayhealing.com

Tel: +86 755-29108212

Fax: +86 755-29108212

Mobile:+86 19820224486

Add: 1st Floor, Building 2, No. 9, Lane 1, Yuanling Industrial Zone,Shangwu

Community, Shiyan Street, Bao'an District, Shenzhen