


# User Manual


---



Model No: CL-M001


# Parts list:

A




x1

B




x1

C




x1

D




x1

E




x1

F




x1

# Product instructions





Tear off the protective film of 3M glue








Attach to proper position as show

# Product instructions



# Touch key description



- 1, To use the special adapter (5V,3A), connect the power supply, press  five seconds to start, you can hear the sound of the buzzer, while the indicator light is on.
- 2, Touch the light respectively  , the buzzer sound is open at the same time, press  and then start running of the product.

## Possible scenarios:

1, When you touch the button, the light will light up and soon go out, this situation is generally.

A: The wire does not match and can not meet the product power output, please change the appropriate power cord to meet the 5V, 3A current;

B: The power output of the adapter is not enough, and it is also necessary to change the adapter to meet the 5V, 3A current.

2, Please do not try to disassemble the product, it will be damaged.

## Contact us

Website: [www.errayhealing.com](http://www.errayhealing.com)

E-mail: [info@errayhealing.com](mailto:info@errayhealing.com)

Tel: +86 755-29108212

Fax: +86 755-29108212

Mobile: +86 19820224486

Add: 1st Floor, Building 2, No. 9, Lane 1, Yuanling Industrial Zone, Shangwu Community, Shiyan Street, Bao'an District, Shenzhen