



Cotton Swabsticks

Pinmed supply Cotton Swabsticks with high quality. Cotton Swabsticks can be used in home and hospital.
*Not for European market

Product Description

Customized Cotton Swabsticks Factory

1. Usage of Cotton Swabsticks

The Cotton Swabsticks is used For Disinfection, topical antiseptic use.

2. Warning of Cotton Swabsticks

- Keep out of reach of children.
- For external use only.
- Do not get into eyes.
- Discard after single use.
- If irritation, please stop use and consult a healthcare provider.
- If accidentally swallowed, seek medical help or contact a Poison Control Center right away.

3. Deliver, Shipping And Serving of Cotton Swabsticks

- 1. Fast Delivery time, (in10-30 days), sometime we have ex-stock too.
- 2. Competitive prices
- 3. OEM&ODM Services(our designer is very professional, we could make design as per your idea for products and packages)
- 4. Full completed documents and certificates for registration.
- 5. Fast communication.
- 6. Quality inspection. Products will be checked during mass production, before out of factory and our QC will check the loading container also.







4.FAQ of Cotton Swabsticks

Q:Are you manufacturer or trading company? A:Industry and trade integration enterprise.

Q:Can i get samples before my order?
A:Free samples are available, freight charge is to customer's account.

Q:Can I get a lower price if I order large quantities?
A:Yes, prices can be discounted with bigger order quantities.

Q:What's the transportation way?
A:DHL,TNT,FEDEX,UPS,EMS, By sea or By Air.

Q:How long it will be to get samples?
A:7-10 days for general products, 15-25days for customized products.

Q:What is the delivery time if I place my order?





A:15-20 days for general products, 25-30 days forcustomized products.

Q:What is the MOQ?

A:No special requirements for MOQ, this can be negotiable.

Q:If OEM is acceptable?

A:Yes,we will try our best to satisfy your demand.